

RBRC 06146

Strain name FVB.Cg-Tg(HSA-Slk*K63R)3405

Primer	primer_name	primer_seq	length	
1		CAGGGCGCTCAGAGTGGACATG	22	RIKEN BRC
2	SLK genoR	CTCTCAGGCGTTAGTGTGCTCTT	24	
3	mLC3ex4AG	GTTAGCATTGAGCTGCAAGCGCCGTCT	27	
4	mLC3ex3GT	TGAGCGAGCTCATCAAGATAATCAGGT	27	


Reaction Components

Enzyme TOYOBO KOD FX (Code.No : KFX-101)

Reaction Conditions

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		0.6			
PCR buffer	2 X	10.0			
dNTP	2.0 mM	4.0			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
primer 4	10 pmol/ μ l	1.0			
DNA		1.0			
Taq	1 U/ μ l	0.4			
total		20			

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	60	30*	
4	68	300*	*Repeat for 30 cycles



Marker

1 kb DNA ladder (Invitrogen, Cat. 15615-024)

Gel

1% agarose in 1 X TAE (100 V, 30 min)

PCR products

primer_set	product_size	product
primer 1 primer 2	ca 3.5 kb	Tg
primer 3 primer 4	550 bp	Internal control (LC3)

Comment Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.

RBRC 06146

Strain name FVB.Cg-Tg(HSA-Slk*K63R)3405

Primer	primer_name	primer_seq	length	
1	TG06-06(HSA)	TGGAATGCGCACAAATATTGAGG	23	RIKEN BRC
2	SLK genoR	CTCTCAGGCGTTAGTGTGCTCTT	24	Depositor's primer
3	mLC3ex4AG	GTTAGCATTGAGCTGCAAGCGCCGTCT	27	Mol. Biol. Cell., 15, 1101-1111(2004).
4	mLC3ex3GT	TGAGCGAGCTCATCAAGATAATCAGGT	27	Mol. Biol. Cell., 15, 1101-1111(2004).


Reaction Components

Enzyme TOYOBO KOD FX (Code.No : KFX-101)

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		0.6			
PCR buffer	2 X	10.0			
dNTP	2.0 mM	4.0			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
primer 4	10 pmol/ μ l	1.0			
DNA		1.0			
Taq	1 U/ μ l	0.4			
total		20			

Reaction Conditions

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	60	30*	
4	68	300*	*Repeat for 30 cycles


Marker

1 kb DNA ladder (Invitrogen, Cat. 15615-024)

Gel

1% agarose in 1 X TAE (100 V, 30 min)

PCR products

primer_set	product_size	product
primer 1 primer 2	ca. 5.0 kbp	Tg
primer 3 primer 4	550 bp	Internal control (LC3)

Comment Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.