



Sheet ID PS_06113
Release date : 12/6/2024
Print date : 12/6/2024


Strain name C57BL/6-Pemt<tm1>

[illegible]

Enzyme	TOYOBO KOD FX (Code.No : KFX-101)
---------------	-----------------------------------

[illegible][illegible]

Gel	2% agarose in 1 X TAE (100 V, 30 min)
------------	---------------------------------------



primer_set		product_size	product (allele type)
primer 1	primer 2	271 bp	WT
primer 1	primer 3	1,094 bp	KO

Comment Standard PCR

Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.