

RBRC 00109

Strain name C3.129S4(B6)-Zic1<tm1Jar>/Jarurbc

[illegible]

Reaction Components


Enzyme TOYOBO KOD FX (Code.No : KFX-101)

Reaction Conditions

[illegible][illegible]

Marker Gene Ladder 100 (Wako, Cat. 316-06951)

Gel	2% agarose in 1 X TAE (100 V, 30 min)
------------	---------------------------------------



primer_set		product_size	product (allele type)
primer 1	primer 2	418 bp	WT
primer 1	primer 3	330 bp	KO

Comment Standard PCR

Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.