

RBRC 00872

Strain name STOCK XO Eda<Ta>/O/Rbrc

Primer	primer_name	primer_seq	length	Reference
1	exon 1	CATGGGCTACCCAGAGGTAGAGC	23	Gene, 371, 42-51(2006).
2	XO_del	CATCTACAAGAGTGTCCAGAACC	25	RIKEN BRC
3	1R2	GGTTGCCCTGTTCCAGTCTGTTTG	25	RIKEN BRC

Reaction Components


Enzyme TOYOBO KOD FX (Code.No : KFX-101)

Reaction Conditions

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	60	30*	
4	68	90*	*Repeat for 30 cycles

Reagent	Conc	Reaction 1	2	3	4
H2O		1.8			
PCR buffer	2 X	10.0			
dNTP	2.0 mM	4.0			
primer 1	10 pmol/μl	1.0			
primer 2	10 pmol/μl	1.0			
primer 3	10 pmol/μl	1.0			
DNA		1.0			
Taq	1 U/μl	0.2			
total		20			

Marker 2-Log DNA Ladder (NEB, Cat. N3200)
Gel 1% agarose in 1 X TAE (100 V, 30 min)



primer_set	product_size	product (allele type)
primer 1	primer 2	1,122 bp WT
primer 2	primer 3	853 bp mutant

Comment Standard PCR

Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.