

RBRC 04200

Strain name B6.129P2-Cxcl12<tm2Tng>

Primer	primer_name	primer_seq	length	
1	Ex1_Fw	GCCGCACTTTCACTCTCGGTCCAC	24	RIKEN BRC
2	intron1_Rv	TCCTATGGCCGCACTTCGCAGAG	23	RIKEN BRC
3	819	CTTGTACAGCTCGTCCATGCCGAG	24	RIKEN BRC

Reaction Components


Enzyme TaKaRa LA Taq with GC buffer (Cat. No. RR02AG)

Reaction Conditions

step	temp (°C)	time (sec)	note
1	94	180	
2	94	30*	
3	60	30*	
4	72	60*	*Repeat for 30 cycles
5	72	300	

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		4.3			
PCR buffer	2 X GC I	10.0			
dNTP	2.5 mM	1.6			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
DNA		1.0			
Taq	5U/ μ l	0.1			
total		20			

Standard PCR



Marker
Gene Ladder 100 (Wako, Cat. 316-06951)

Gel
2% agarose in 1 X TAE (100 V, 30 min)

PCR products			
	primer_set	product_size	product
primer 1	primer 2	379 bp	WT
primer 1	primer 3	ca 900 bp	KI

Comment