

RBRC 05691

Strain name B6;129-Fgf18<tm1Tri>

Primer	primer_name	primer_seq	length	
1	int3f	CTGAGGAGTGTGGACAGGGAGG	22	RIKEN BRC
2	ex3rv	CATTCTCCTCGGCTGCCAACAC	22	RIKEN BRC

**Reaction Components**


Enzyme TOYOBO KOD FX (Code.No : KFX-101)

**Reaction Conditions**

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	68	120*	*Repeat for 30 cycles

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		2.8			
PCR buffer	2 X	10.0			
dNTP	2.0 mM	4.0			
primer 1	10 pmol/ $\mu$ l	1.0			
primer 2	10 pmol/ $\mu$ l	1.0			
DNA		1.0			
Taq	1 U/ $\mu$ l	0.2			
<b>total</b>		<b>20</b>			

Standard PCR



**Marker**  
Gene Ladder 100 (Wako, Cat. 316-06951)

**Gel**  
2% agarose in 1 X TAE (100 V, 30 min)

**PCR products**

primer_set	product_size	product
primer 1 primer 2	704 bp	WT
primer 1 primer 2	ca 800 bp	flox

**Comment** Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.