

RBRC 06156

Strain name C57BL/6N-Slc5a4a<tm1Csk>

Primer	primer_name	primer_seq	length	Reference
1	SGLT3a wild F2	TGAGGCTGATCATTAACTCAGGAC	24	Depositor's primer
2	SGLT3a wild R2	TTCGGTTCGTCTTCAGCATTGC	22	Depositor's primer
3	Neo-anti4	AACTTCCTGACTAGGGGAGGAGTA	24	Depositor's primer

Reaction Components


Enzyme TOYOBO KOD FX (Code.No : KFX-101)

Reaction Conditions

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	60	30*	
4	68	60*	*Repeat for 30 cycles

Reagent	Conc	Reaction 1	2	3	4	(μL)
H2O		1.8				
PCR buffer	2 X	10.0				
dNTP	2.0 mM	4.0				
primer 1	10 pmol/μl	1.0				
primer 2	10 pmol/μl	1.0				
primer 3	10 pmol/μl	1.0				
DNA		1.0				
Taq	1 U/μl	0.2				
total		20				

Marker Gene Ladder 100 (Wako, Cat. 316-06951)
Gel 2% agarose in 1 X TAE (100 V, 30 min)



primer_set	product_size	product (allele type)	
primer 1	932 bp	WT	+/+ : 932 bp +/- : 932 bp, 256 bp
primer 1	256 bp	KO	

Comment Standard PCR

Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.