

RBRC 06330

Strain name C57BL/6-Tg(LY6A-KDM1A\*v2)6Fscr

Primer	primer_name	primer_seq	length	
1		TGGCTTGATCAGGTCACAAA	20	Depositor's primer
2		TCTCTTAGGAACAGCTTGTCCG	23	RIKEN BRC
3	mLC3ex4AG	GTTAGCATTGAGCTGCAAGCGCCGTCT	27	Mol. Biol. Cell., 15, 1101-1111(2004).
4	mLC3ex3GT	TGAGCGAGCTCATCAAGATAATCAGGT	27	Mol. Biol. Cell., 15, 1101-1111(2004).


### Reaction Components

Enzyme TOYOBO KOD FX (Code.No : KFX-101)

### Reaction Conditions

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		0.8			
PCR buffer	2 X	10.0			
dNTP	2.0 mM	4.0			
primer 1	10 pmol/μl	1.0			
primer 2	10 pmol/μl	1.0			
primer 3	10 pmol/μl	1.0			
primer 4	10 pmol/μl	1.0			
DNA		1.0			
Taq	1 U/μl	0.2			
<b>total</b>		<b>20</b>			

step	temp (°C)	time (sec)	note
1	94	120	
2	98	10*	
3	60	30*	
4	68	120*	*Repeat for 30 cycles



### Marker

1 kb DNA ladder (Invitrogen, Cat. 15615-024)

### Gel

1% agarose in 1 X TAE (100 V, 30 min)

### PCR products

primer_set	product_size	product
primer 1 primer 2	1,439 bp	Tg
primer 3 primer 4	550 bp	Internal control (LC3)

**Comment** Novagen KOD Xtreme™ Hot Start DNA Polymerase can be purchased from EMD.